

TorqTech handpieces

Models and areas of use

1:5



TorqTech CA-5IF-O

Use: E.g. for preparation work with 1.6 mm FG shaft type diamond dental burs

1:1



TorqTech CA-DC-O

Use: E.g. for working of dentine and dentinal matter with 2.35 mm RA shaft type carbide-metal instruments

1:1



TorqTech ST-DH

Use: E.g. for corrective work on dental prostheses with 2.35 mm HP shaft type grinding instruments

4:1



TorqTech CA-4R-PP

Use: E.g. for polishing with screw-in rubber wedge or small brush

TorqTech handpieces

Models and areas of use

4:1



TorqTech CA-4R-O-PF

Use: E.g. for polishing with EVA® Tips or for fine grinding with Lamineer® Tips

10:1



TorqTech CA-10RC-O

Use: E.g. for surgical work with 2.35 mm RA shaft type carbide-metal instruments

10:1



TorqTech CA-10RC

Use: E.g. for mechanical root canal preparation work using nickel-titanium instruments

10:1



TorqTech CA-10RC-Endo*

Use: For safe root canal preparation work using nickel-titanium instruments while applying constant length control

* Endo contra angle only for Soaric