

Press release

## The new Veraviewepocs 3D R100 of Morita wins Innovations Award 2012

Dietzenbach, July 2012. The latest Morita CBCT scanner won the first Innovations Award at the Clinical Innovations Conference in London. The prize showcases the most pioneering product on the dental market. The Veraviewepocs 3D R100 received the award for its revolutionary Reuleaux triangle full arch field of view. By closely matching the natural dental arch, patient x-ray dosage is reduced by around 15 percent, compared with a conventional 100 mm cylindrical scan.

**Contact:**

J. Morita Europe GmbH

Julia Meyn

Justus-von-Liebig-Straße 27a

63128 Dietzenbach

Germany

T +49. 6074. 836 110

F +49. 6074. 836 299

[jmeyn@morita.de](mailto:jmeyn@morita.de)

[www.morita.com/europe](http://www.morita.com/europe)