

Dental Simulator
SSX 2.5




Slide Arm Type

- ▷ With Warning Light
- ▷ With Out Cabinet




Specifications

Model Designation	: SSX 2.5
Input Voltage	: 220V AC $\pm$ 10%
Output Voltage	: 0-30V AC
Frequency	: 50 / 60 Hz
Apparent power	: 1.5 kVA
Current	: 6A at 220V
Main fuse	: 250V AC 8A
Phantom & Torso	: Max 780 $\pm$ 20 mm
Height	: Min 680 $\pm$ 20 mm
Phantom & Torso	: Max 75 $\pm$ 2°
Tilt Angle	: Min 0 $\pm$ 2° (From Horizontal Level)
Weight	: 132 kg.



* Dimension tolerance $\pm$ 20 mm



Standard Equipment & Feature

LED Operating Light

- On/Off By Automatic Positioning
- On/Off By Foot Control Switch
- Brightness Adjust By Panel Switch

Warning Light

- Switch For Warning Light On/Off

Doctor side

- Doctor Tray table with 5 Instrument Holders
- High Speed Air Turbine Circuit With Fiber Optic 4H Tubing (1)
- Low Speed Air Motor Circuit 4H Tubing (1)
- Doctor Side 3-Way Syringe With Tubing (1)

SE Holder & Flex Arm

- Saliva Ejector With Tubing (1)

Phantom & Torso

- Electric Motor Drive
- Manual Operation Control By Panel Switch
- Two preset positioning & Reset position Control By Panel Switch

Foot Controller

- Foot Pedal For Handpieces
- Switch For Light On/Off
- Switch For Water On/Off

Vacuum System

- Air Vacuum System

SIAMDENT.

A Morita Group Company

Siamdent co., Ltd.

71/10 Moo. 5, T. Tharkham, A. Bangpakong, Chachuengsao 24130 Thailand

TEL:+66 (0) 3857 3042, FAX:+66 (0) 3857 3043

www.siamdent.com

S-EDS-024-0017/04-08-26/A